ergositter

The back recovery chair!

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The Jury's Motivations

The winner of the Life Science category: Backrefresh seating system WO- BPS. A simple and daring take on back pain, a major medical problem. The solution designed by the 2016 is a system that is well-substantiated by scientific and medical findings, in which an ingenious chair relieves pressure and elongate the sitter's back to both cure and prevent disc herniation.

Karin Wanngård, Mayor of Stockholm, Sweden

"I assumed that something like this would show up sometime in the future, but not this soon!"

Gunnar Andersson, MD, Professor of Orthopedics, Chicago, USA. One of the world's foremost back specialists and an expert on sitting physiology.

"This is the way all people will be sitting in the future."

Christian Höfors Ziebell, PhD, ASSCO prof. Chalmers. Biomechanic, back and shoulder specialist.

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100% ryggbelastning sittende i vanlig stol.

<0% ryggbelastning sittende i Ergositterposition.

Sitting with Adjustable Ischial and Back supports: Biomechanical changes

Mohsen Makhsous, PhD^{a,b}, Fang Lin, PhD^{a,b}, Ronald W. Hendrix, MD^e, Matthew Hepler, MD^c, and Li-Qun Zhang, PhD^{a,d}

^aSensory Motor Performance Program, Rehabilitation Institute of Chicago, ^bDepartments of Physical Medicine & Rehabilitation, ^cOrthopaedic Surgery, ^dBiomedical Engineering, and ^eRadiology, Northwestern University

INTRODUCTION: Sitting may induce posterior rotation of pelvis, reduction of lumbar lordosis, and increases in muscle tension, disc pressure, and pressure on ischium and coccyx, which may be associated with low back pain (LBP) [1-2].

- A device (Fig. 1) that reduces ischial load and maintains lumbar lordosis may help increase seating comfort and reduce LBP.

- Objective: To investigate the biomechanical effects of adjusting ischial and/or backrest supports during sitting.

METHODS: 15 healthy office workers

- Sitting conditions: sitting with & without ischial support & with adjustable back support
- Sacral inclination, Lumbar lordosis, Intervertebral space (Fig. 2)
- Muscular activity in stabilizing the trunk
- Contact pressure distributions (Fig. 3)
- Reaction forces between buttock-thighs & seat between back and backrest
- Load carried by seat pan and backrest


Fig 3. Representative pressure distribution on the seat (the lower row) and backrest (upper row) for 'Upright with backrest' sitting condition.

Left column: 'With Ischial Support' when the BPS was at level = 0 position.

Right column: 'Without Ischial Support' when the BPS was tilted down to = 18 .

RESULTS:

When the ischial support was relieved, the center of the force on the seat and on the legs of the chair, and the peak center of pressure on the seat, were significantly ($P < 0.002$) shifted forward towards the thighs.

- The total contact area on the seat pan and on the backrest was significantly decreased and increased, respectively ($P < 0.001$).

- The sacral inclination, total & segmental lumbar lordosis, and lumbar spine disc height were significantly increased for sitting 'Upright with Backrest' with the lumbar curve close to that during standing. (Table 1)


Fig 1. The seat of the chair consists of two parts, i.e. FPS, which is fixed and BPS, which can be tilted downward or upward by = 18 . The shape of the low back support can be adjusted for individual subject using an inflatable air-filled cushion.


Fig 2. Representative sacral inclination (), total and segmental lumbar lordosis from one subject for sitting 'Upright with fully fitted Backrest' with the BPS at level (= 0) (A) and tilted down (= 18) (B) position, and at the standing (C) postures.

References: [1].Andersson GJB, et al. The influence of backrest inclination and lumbar support on lumbar lordosis. Spine 1979; 4:52-8. [2]. Lord MJ, et al. Lumbar Lordosis. Effects of sitting and standing. SPINE 1997; 22:2571-4.

Table 1. The sacral inclination (), total (S1-L1) and segmental (S1-L2, S1-L3, S1-L4, S1-L5) lumbar lordosis for standing and sitting postures.

| Posture | Sitting 'Upright with Backrest' | | Standing |
|----------------------------|---------------------------------|--------------------|--------------------|
| | Fitted partially (n=6) | Fitted fully (n=9) | |
| Ischial support | With | Without | |
| Sacral inclination: | | | |
| α | 14.04 \pm 11.30° | 20.95 \pm 11.22° | 14.14 \pm 10.94° |
| P | | 0.001 | 24.95 \pm 9.50° |
| | | | <0.001 |
| Total Lordosis: | | | |
| S1-L1 | 24.53 \pm 21.31° | 32.63 \pm 22.25° | 38.94 \pm 14.61° |
| P | | 0.010 | 47.86 \pm 11.89° |
| | | | 0.008 |
| Segmental Lordosis: | | | |
| S1_L2 | 25.26 \pm 16.48° | 33.06 \pm 17.96° | 36.04 \pm 9.24° |
| P | | 0.015 | 42.92 \pm 9.75° |
| | | | 0.005 |
| S1_L3 | 22.10 \pm 12.92° | 30.65 \pm 15.61° | 29.71 \pm 6.14° |
| P | | 0.023 | 35.20 \pm 10.88° |
| | | | 0.041 |
| S1_L4 | 16.31 \pm 8.58° | 24.62 \pm 12.41° | 21.15 \pm 5.56° |
| P | | 0.024 | 27.53 \pm 9.09° |
| | | | 0.025 |
| S1_L5 | 12.41 \pm 5.96° | 19.10 \pm 6.45° | 11.12 \pm 3.88° |
| P | | 0.006 | 17.74 \pm 4.79° |
| | | | 0.001 |

SUMMARY & CONCLUSIONS: Sitting with reduced ischial support and fitted backrest to the lower spine

- altered the contact area
 - resulted in reduced peak pressure under the ischia
 - reduced muscular activity
 - maintained total and segmental lumbar lordosis
 - forwardly rotated the sacrum
 - increased lumbar intervertebral disc heights
- which could potentially reduce LBP.

ACKNOWLEDGEMENT: The project was supported in part by NIDRR (LZ) and Falk Medical Research Trust (LZ)

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